




3M™ Health Care Respirators



EN 149 Classification	Catalogue No.	Valved/ Unvalved	Minimum % Total Efficiency	Maximum Inward Leakage %	
FFP1S	1861	Unvalved	78	22	
FFP2S	1862	Unvalved	92	8	
FFP2S	1872V	Valved	92	8	
FFP3S	1863	Unvalved	98	2	
FFP3S	1873V	Valved	98	2	

- Protects against airborne micro-organisms, laser plume and diathermy smoke
- Complies with EC Personal Protective Equipment Directive 89/686
- EN149 approved - FFP1S, FFP2S, or FFP3S
- Low profile, three panel design
- Valved or unvalved

3M Health Care

3M is a trademark of the 3M Company.
© 3M Health Care 2003.

3M Health Care Limited

3M House
Morley Street
Loughborough
Leicestershire LE11 1EP
Tele: (01509) 613151
Website: www.3mhealthcare.co.uk



3M™ Infection Prevention